# BT110

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# BT110

## Self-Contained Biological Indicator.

*Rapid Readout Fluorescence System.*

### Usage

Monitoring Ethylene Oxide Sterilization Processes.

### Applicable Regulation

ISO 11138-1:2006 and ISO 11138-2:2006; IRAM 37102-1:1999 and IRAM 37102-2:1999.

### Classification

Class 1, according to risk.

### Authorization

Designed under Quality Management System standards ISO 9001:2008 and ISO 13485:2012. ANMAT (Argentinean National Administration of Drugs, Food and Medical Technology) PM 1614-1.

### Characteristics

Polypropylene random copolymer tube: 50.0 mm high x 8.5 mm external diameter. Wall thickness: 0.5 mm

Polypropylene microfiber cap at the bottom of the tube.

Polypropylene cap: 10.5 mm external diameter, 16.5 mm high. Wall thickness: 1.0 mm

Glass ampoule: 35.0 – 38.0 mm high. mm. External diameter: 6.5 mm. Wall thickness: 0.2 – 0.3 mm.

Culture medium 0.5 – 0.7 ml, **blue** color.

Cap filter: medical grade paper, 60 g/m<sup>2</sup>, 16.0 mm diameter.

Filter Paper: 5.0 x 20.0 mm, 160 g/m<sup>2</sup>.

≥ 10<sup>6</sup> *Bacillus atropheus* ATCC 9372 spores per vial.

Final fluorescence reading: at 4 hour-incubation at 37 °C (sensitivity: 97%).

An optional visual pH color change confirmation could be made after 48 hours of incubation. If sterilization process has not been successful, culture medium will change to a greenish color first, and then to yellow during incubation at 37 °C, thus showing the presence of living spores. If sterilization process is successful culture medium will remain blue after the incubation process.

7-day readout for visual confirmation by color change is optional and not intended to be routinely performed; it is an initial validation of the 4 hour-reading. Fluorescence results may be compared to the 7-day visual reading. If sterilization process has not been successful, culture medium will change to a greenish color and then to **yellow** during incubation at 37 °C, thus showing the presence of living spores. If sterilization process is successful culture medium will remain blue after the incubation process.

NOTE: if 7-day readout is performed, a humidified environment will be required to avoid medium dry out.

D-Value: Not lower than 2.5 minutes, conditions: 600 ± 30 mg/l EO, 60 ± 10 % RH, 54 ± 1 °C.

### Environmental conditions during manufacture

Room temperature 15-30 °C, RH 30-80 %. Sterility conditions are necessary only during the inoculation process performed in laminar flow.

### Storage conditions

Best stored in original box under the next conditions: T = 10-30 °C, RH 30-80 %.

### **Transportation conditions**

Storage conditions should be strictly followed.

Products should be transported in closed and reinforced boxes in order to avoid damages. Product transportation does not represent any risk for human health.

### **Shelf-life**

24 months.

### **Packing**

50 units per box.

Packing information: Product code and description, process for intended use, presentation, classification according to regulation, manufacturer information and data on box's label. Note: manufacture date is calculated by subtracting 24 months to the expiration date.

Weight per box: 143,7 g.

### **Labelling**

On product: 17.0 mm x 33.0 mm polypropylene label. Printed in black. 1.5 mm chemical indicator line, printed with EO reactive ink (Color change to green). Product code and batch number, process for intended use and bacterial strain.

On product's box: product code, batch number, bacterial load, manufacture and expiration date.

### **Possible target markets**

Healthcare and Industry.

### **Other important information (optional)**

It is advisable to incubate at 37 °C in Bionova® IC10/20FR Auto-Reader incubator.

Read product's directions for use thoroughly before use.

**Note:** When necessary, limits of each product's parameters should be specified.

### **Description of identified and/or predictable risks which could represent a risk for the product**

Do not store the product near sterilizing agents.

Do not expose this product to Steam, Dry Heat, Radiation or any sterilization process other than EO.

The positive biological indicator can be autoclaved at 121 °C for at least 20 minutes, at 132 °C for at least 15 minutes in gravity displacement steam sterilization, or at 134 °C for at least 10 minutes in vacuum steam sterilizer.