IT26-C

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IT26-C

Moving Front Integrator Indicator

For Steam sterilization processes.

Usage

Control of Steam sterilization processes between 118 °C and 138 °C.

Applicable Regulation

EN ISO 11140-1:2009

Classification

Class 1, according to risk.

Authorization

Designed under Quality Management Systems ISO 9001:2008 and ISO 13485:2012. ANMAT (Argentinean National Administration of Drugs, Food and Medical Technology) PM 1614-4.

Characteristics

Class 5 Integrator Indicator.

Integrator Indicator consists of a paper wick and a steam and temperature sensitive chemical pellet contained in a paper/film/foil laminate. The chemical pellet melts and migrates as a dark bar along the paper wick. The migration is visible through a zone marked ACCEPT or REJECT, thus indicating whether sterilization conditions were met. The extent of migration depends on steam quality, time, and temperature. The ACCEPT result is reached when a theoretical spore population reaches its kill time, indicating integration condition has been reached. This condition is calibrated with the kill time of a *10⁶Geobacillusstearotherophilus* spores population, calculated in BIER (Biological Indicator Evaluator Resistometer).

Size: 52 mm x 21 mm.

Conditions: saturated steam at 121 °C, 128 °C, 135 °C.

Environmental conditions during manufacture

15-30 °C (room temperature), RH 30-80 %.

Storage conditions

T = 10-30 °C, RH 35-70 %, keep out of direct light.

Transportation conditions

Storage conditions should be strictly followed. Transport in closed and reinforced boxes in order to avoid damages. Transportation of this product does not represent a risk for health.

Shelf-life

5 years

Packing

100 units per bag.

Packing information: product code and description, process for intended use, operating instructions, product presentation, classification according to regulation, batch number, manufacturer information, expiration date and SV at 121 °C, 128 °C and 135 °C.

Weight per package: 88,1 g

Labelling

Information only on product's box.

Possible target markets

Healthcare, Food, Pharmaceutical and Medical Industries.

Other important information (optional)

Read product's directions for use thoroughly before use.

Note: When necessary, limits for each product's parameters should be specified.

Description of identified and/or predictable risks which could represent a risk for the product

Do not store the product near sterilizing agents.

Do not expose this product to EO, Dry Heat or any sterilization processes other than Steam.

Do not reuse the sterilizer until the dark bar has reached the ACCEPT zone.

Do not store near strong alkaline or acidic products such as cleaning/disinfecting agents.