CDWU

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CDWU

Ultrasonic Cleaning Indicator

For ultrasonic cavitation performance tests.

Usage

Designed for ultrasonic cleaning process routine control and monitoring ultrasonic cavitation effectiveness of ultrasonic washing machines operating at frequencies equal or greater than 35 kHz, at temperatures ranging from 18 °C to 70 °C.

Authorization

Designed under Quality Controls integrated systems ISO 9001:2008 and ISO 13485:2012.

Characteristics

2 ml clear glass (borosilicate) vial: 31.0 mm high x 11.6 mm external diameter, 8.0 mm neck diameter. Wall thickness: 0.5 mm.

Brown/Red viton septa: 8.0 x 1 mm.

Aluminum top cap: 5.45 mm high x 11.66 mm diameter. Wall thickness: 0.34 mm.

Glass (borosilicate) beads: 3.0 mm.

Blue clear reactive solution: 1.5 ml.

Allows detection of deficiencies in the cavitation performance of the ultrasonic washing machine.

Challenge level: high

Initial Color: blue

Result: yellow transparent solution

NOTE: If unfavorable results are obtained, check working frequency of the ultrasonic cleaner, time and temperature of the ultrasonic washing, positioning and amount of wash load, presence of bubbles within the wash solution, location of indicator and general condition of ultrasonic washing machine.

Environmental conditions during manufacture

10-30 °C (room temperature), RH 30-80%.

Storage conditions

T = 10-30 °C, RH 30-80%, store in dark place.

Transportation conditions

Storage conditions should be strictly followed.

Transport in closed and strengthened boxes in order to avoid damages. Transportation of this product does not represent a risk for health. Product can be kept at temperature below 25 °C for up to 15 days during transportation.

Shelf-life

1 year

Packing

30 units per box.

Packing information: product code and description, process for intended use, instructions for use, manufacturer information and data on label.

Labelling

On product's box: Product code and description, process for intended use, batch number and expiration date.

Possible target markets

Healthcare and Industry.

Other important information (optional)

Read product's directions for use thoroughly before use.

NOTE: When necessary, limits for each product's parameters should be specified.

Description of identified and/or predictable risks which could represent a risk for the product

CDWU indicators should be stored at temperatures between 10-30 °C, specially protected from light and moisture, within its original packaging.

Do not re-use the CDWU indicator.

Do not open the indicator before its use, since this will cause its deterioration.

Do not use the CDWU indicator for monitoring normal washing using non-ultrasonic cleaners.