KPCD222-C

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KPCD222-C

Steam Test Pack PCD Kit

Process Challenge Device for Steam sterilization processes

Usage

For routine load release, specially for loads containing implants. For routine sterilizer monitoring. For periodic validation of the sterilizer (after installation, repair, relocation)

Detects inadequate air removal and steam penetration in dynamic-air-removal (pre-vacuum) and gravity displacement steam autoclaves.

Condition: 132/135 °C ≥ 4 minutes and 121 °C ≥ 30 minutes, steam sterilization processes.

Applicable Regulation

EN ISO 11140-1:2009, ISO 11138-1:2006, ISO 11138-3:2006, ANSI/AAMI ST79

Classification

Class 1, according to risk.

Authorization

Designed under Quality Management Systems ISO 9001:2008 and ISO 13485:2012. ANMAT (Argentinean National Administration of Drugs, Food and Medical Technology) PM 1614-4.

Characteristics

Porous cards system holding a Self-contained Biological Indicator and IT26-C moving-front Class 5 Integrator Indicator. The system also contains a self-adhesive Record Card, PCDBI-C-RC, where data on sterilization cycle may be written.

Cardboard box (10.5 cm x 6.6 cm x 1.8 cm) with a self-adhesive Class 1 Process Indicator.

Process Indicator:

Initial color: **blue**

Final color: **dark grey or black**

IT26-C Integrator Indicator (within the challenge device):

The migration is visible through a zone marked ACCEPT or REJECT, thus indicating whether sterilization conditions were met.

Integration condition is calibrated against the kill time of a 10^6 *G. stearothermophilus* ATCC 7953 spores' population, calculated in BIER (Biological Indicator Evaluator). Conditions: saturated steam at 121 °C, 128 °C, 135 °C.

Self-contained Super Rapid Biological Indicator (within the challenge device): ≥ 10^6 *G. stearothermophilus* ATCC 7953 spores per vial in paper carrier. Glass ampoule with growth indicator medium. Label with chemical indicator line printed with Steam reactive ink (Color change: pink to brown).

Final results: by fluorescence after incubation for 1 hour at 60 °C (sensitivity 97%). A 7-day readout is optional to confirm 1-hour results by color change.

NOTE: If sterilization process was not successful, indicator medium will turn yellow after incubation, thus indicating the presence of living *G. stearothermophilus* spores. If sterilization process was successful, indicator medium will remain purple after incubation.

Environmental conditions during manufacture

Room temperature 15-30 °C, RH 30-80 %.

Storage conditions

T = 10-30 °C, RH 35-70 %, keep in a dark place.

Transportation conditions

Storage conditions should be strictly followed. Carry in closed and reinforced boxes in order to avoid damages. Transportation of this product does not represent a risk for health.

Shelf-life

2 years

Packing

25 PCD222-C + 25 BT222 Biological Indicators per box.

Packing information: product description and manufacturer information.

Labelling

On product's box: product code and description, process for intended use, batch number and expiration date.

Possible target markets

Healthcare, Food, Pharmaceutical and Medical Industries.

Other important information (optional)

Read product's directions for use thoroughly before use.

Note: When necessary, limits for each product's parameters should be specified.

Description of identified and/or predictable risks which could represent a risk for the product

Do not store the product near sterilizing agents.

Do not expose this product to EO, Dry Heat, Radiation or any sterilization process other than Steam.

The positive biological indicator can be autoclaved at 121 °C for at least 20 minutes, at 132 °C for at least 15 minutes in gravity displacement steam sterilization, or at 134 °C for at least 10 minutes in vacuum steam sterilizer.