BT220

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BT220

Self-Contained Biological Indicator.

Rapid Readout Fluorescence System.

Usage

Monitoring vacuum-assisted Steam Sterilization processes at 121 °C – 135 °C and gravity air-displacement Steam Sterilization processes at 121 °C – 135 °C.

Applicable Regulation

ISO 11138-1:2006 and ISO 11138-3:2006; IRAM 37102-1:1999 and IRAM 37102-3:1999.

Classification

Class 1, according to risk.

Authorization

Designed under Quality Management System standards ISO 9001:2008 and ISO 13485:2012. ANMAT (Argentinean National Administration of Drugs, Food and Medical Technology) PM 1614-1.

Characteristics

Polypropylene random copolymer tube: 50.0 mm high x 8.5 mm external diameter. Wall thickness: 0.5 mm

Polypropylene microfiber cap at the bottom of the tube.

Polypropylene cap: 10.5 mm external diameter, 16.5 mm high. Wall thickness: 1.0 mm

Glass ampoule: 35.0 – 38.0 mm high. mm. External diameter: 6.5 mm. Wall thickness: 0.2 – 0.3 mm.

Culture medium 0.5 – 0.7 ml, purple color.

Cap filter: medical grade paper, 60 g/m², 16.0 mm diameter.

Spores carrier: Filter Paper: 5.0 x 20.0 mm, 160 g/m².

≥ 10⁶ *Geobacillus stearothermophilus* ATCC 7953 spores per vial.

Final fluorescence reading is performed after 3 hour-incubation at 60 °C (sensitivity: 97 %).

An optional visual pH color change confirmation could be made after 48 hours of incubation. If sterilization process has not been successful, culture medium will change to a greenish color first, and then to yellow during incubation at 60 °C, thus showing the presence of living spores. If sterilization process is successful culture medium will remain purple after the incubation process.

7-day readout is optional and not intended to be routinely performed; it is an initial validation of the 3 hour-reading. Fluorescence results may be compared to the 7-day visual reading.

NOTE: If 7-day readout is performed, a humidified environment will be required to avoid medium to dry out.

D-Value: Not lower than 1.5 minutes at 121 °C. Another D-value is informed at 135°C.

Z-Value: Not lower than 6 °C.

Environmental conditions during manufacture

Room temperature 15-30 °C, RH 30-80 %. Sterility conditions are necessary only during inoculation process performed in laminar flow.

Storage conditions

Best stored in original box under the next conditions: T = 10-30 °C, RH 30-80 %.

Transportation conditions

Storage conditions should be strictly followed.

Products should be transported in closed and reinforced boxes in order to avoid damages. Product transportation does not represent any risk for human health.

Shelf-life

24 months.

Packing

50 units per box.

Packing information: product code and description, process for intended use, presentation, classification according to regulation, manufacturer information and data on box's label.

Weight per box: 142,2 g.

NOTE: manufacture date is calculated by subtracting 24 months to the expiration date.

Labelling

On product: 17.0 mm x 33.0 mm polypropylene label. Printed in black. 1.5 mm chemical indicator line, printed with Steam reactive ink (Color change to brown). Product code and batch number, process for intended use and bacterial strain.

On product's box: product code, batch number, bacterial load, manufacture and expiration date.

Possible target markets

Healthcare and Industry.

Other important information (optional)

It is advisable to incubate at 60 °C in Bionova® IC10/20FR Auto-Reader incubators.

Read product's directions for use thoroughly before use.

NOTE: When necessary, limits of each product's parameters should be specified.

Description of identified and/or predictable risks which could represent a risk for the product

Do not store the product near sterilizing agents.

Do not expose this product to EO, Dry Heat, Radiation or any sterilization process other than Steam.

The positive biological indicators can be autoclaved at 121 °C for at least 20 minutes, at 132 °C for at least 15 minutes in gravity displacement steam sterilization, or at 134 °C for at least 10 minutes in vacuum steam sterilizer.