



NEWSTERI SRL

Sterilization products.

EtoPro[®]
ethylene oxide specialists

Product Brochure

Cartridge



ETOPRO ® offers aluminium cartridges of the highest purity in controlled quantities and packaged under severe safety regulations which fulfill all the requirements of the sterilizers present in the whole market.

- One use cartridges
- High effectiveness as a bactericide, fungicide and antivirus
- Compatible with any sterilizer of the market

ETOPRO cartridges are produced 100% with European supplies ensuring high quality and a life of two years

Information

EtoPro ® cartridges				
Code	Sizes		Ethylene Oxide content	
	Diameter	Height	Gram	Ounces
EP-25	25 mm	105 mm	30 grs	1,07 oz
EP-50	35 mm	105 mm	67 grs	2,37 oz
EP-67	35 mm	105 mm	50 grs	1,77 oz
EP-70	35 mm	140 mm	70 grs	2,48 oz
EP-100	38 mm	165 mm	100 grs	3,52 oz
EP-120	50 mm	133 mm	120 grs	4,24 oz
EP-130	50 mm	133 mm	130 grs	4,60 oz
EP-134	38 mm	165 mm	127 grs	4,47 oz
EP-170	45 mm	178 mm	170 grs	5,99 oz
EP-200	50 mm	174 mm	200 grs	7,07 oz
EP-230	50 mm	174 mm	230 grs	8,14 oz

Ampoules



Sterilization through Ethylene oxide ampoules ETOPRO is a simple method which can be practiced anywhere specially in critical conditions. This method allows to sterilize any type of material including those which can not be exposed to sterilization with steam or dry heat.

This method assures that the user will achieve a correct sterilization process.

The ETOPRO kit has all the necessary elements to acquire an appropriate sterilization.

What it contains?

- ▶ Ethylene Oxide ampoule in a recipient with a diffusion membrane
 - ▶ Humidigene
 - ▶ Integraton class 5, Integron
- ▶ Contention bag of the material which will be exposed to sterilization
 - ▶ Seal to close the bag

Information

Codes	Volume	ETO	Bag	Capacity	Box
EP-5	5grs	90 %	31 x 60 cm	11 lts	60 pieces
EP-20	20grs	97 %	55 x 90 cm	34 lts	15 pieces